



Pocket magnifier

ISOLAB

Pocket magnifier

- Suitable for reading small prints with a magnification of 2X.
- 3 LED bulbs with on/off control.
- Ergonomic design
- Protective bag included



Code	Description	Packaging
LLG06267408	Pocket magnifier, 2X magnification	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

