



## Compressed air powered stirrers, PLR, accessories

### Compressed air powered stirrers, PLR, accessories

#### Connection couplings VK:

Fixed-diameter connection for motor/magnetic coupling drive shaft and stirring rod. Suitable for use in potentially explosive atmospheres. Material No. 1.4404.

#### Flexible glass stirring shaft couplings PK:

For stirring tasks using glass stirring rods. Elastic connector between motor and stirring shaft. Provides exact alignment of motor and shaft, and evens out axial and height differences. Accepts 6mm to 10mm diameter shafts.

#### Flexible glass stirring shaft couplings K:

For stirring tasks using glass stirring rods. Elastic connector between motor and stirring shaft. Provides exact alignment of motor and shaft, and evens out axial and height differences. Accepts 6mm to 10mm diameter shafts

#### Double-Cardan-Coupling DCK I:

With 8mm dia. hexagonal input drive receiver and 6mm square socket drive output.



Code	Description	Packaging
LLG09779052	Compressed air powered stirrers, PLR, accessories, Type F 8, Chuck for shafts 0.5 to 8.5mm	1 pz.
LLG06285834	Compressed air powered stirrers, PLR, accessories, Type VK 10x10, Connection couplings 10 x 10 mm	1 pz.
LLG06286426	Compressed air powered stirrers, PLR, accessories, Type VK 10x8, Connection couplings 10 x 8 mm	1 pz.
LLG09197400	Compressed air powered stirrers, PLR, accessories, Type K 10/13, Flex coupling receiver i. d. 10/13mm	1 pz.
LLG09779065	Compressed air powered stirrers, PLR, accessories, Type DCK I, Stainless steel, L = 110mm	1 pz.
LLG09198010	Compressed air powered stirrers, PLR, accessories, Type K 10, Flex coupling receiver i. d. 10mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

