


**Multi-Rack for syringe pumps Legato 200 and Legato 210****KD Scientific****Multi-Rack for syringe pumps Legato® 200 and Legato® 210**

Modular syringe racks can be purchased to create a multichannel syringe pump. These racks will work with the Legato® 200 or Legato® 210 only. All of them could be for infuse only as well.

Small Syringe Multi-Rack: The small syringe multi-rack will accommodate up to six 30ml to 60ml syringes or up to ten 0.5µl to 20ml syringes. For Infuse/Withdraw applications.

Large Syringe Multi-Rack: The large syringe multi-rack will accommodate up to four 60ml to 140ml plastic syringes. For Infuse/Withdraw applications.

Microliter Syringe Multi-Rack: The microliter syringe multi-rack will accommodate up to four 0.5µl to 10ml syringes. For Infuse/Withdraw applications.



Code	Description	Packaging
LLG06267704	Multi-Rack for syringe pumps Legato® 200 and Legato® 210, Large Syringe Multi-Rack	1 pz.
LLG06267705	Multi-Rack for syringe pumps Legato® 200 and Legato® 210, Microliter Syringe Multi-Rack	1 pz.
LLG06267703	Multi-Rack for syringe pumps Legato® 200 and Legato® 210, Small Syringe Multi-Rack	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

