

**Rack for centrifuge tubes, PP****ISOLAB****Rack for centrifuge tubes, PP**

Suitable to hold 15 ml and 50 ml centrifuge tubes with conical bottom. Just push the tubes into the cavities and leave them. Unique design has an angular end to provide easy stacking for space saving purpose. Moulded in alpha numeric grid reference on the rack allows easy sample identification.



Code	Description	Packaging
LLG06267412	Rack for centrifuge tubes, PP, For tubes 15 ml, Array 5 x 5, dim. (W x D x H) 150 x 120 x 32 mm	1 pz.
LLG06267413	Rack for centrifuge tubes, PP, For tubes 50 ml, Array 5 x 5, dim. (W x D x H) 230 x 190 x 35 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

