

**LLG-Cryotubes, PP, sterile****LLG Labware****LLG-Cryotubes, PP, sterile**

With star-shape in the screw cap, for automated capper/decapper.

Temperature resistant down to  $-196^{\circ}\text{C}$ . **Not appropriate for operations in liquid nitrogen.** With marking area, fill line and graduations printed on the vial. Safe and precise handling of biotechnology materials. Bar-code printed on each tube. DNase/RNase free. Vials and caps are autoclavable at  $121^{\circ}\text{C}$ .

**Homogeneous conception:**

-All made from polypropylene

-Single-body, single-turn screw cap

Printed graduations for accurate measurements. Large white area for writing specimen identification. Sterile.

Packed in 2 bags of 50 tubes.



Code	Description	Packaging
LLG09401161	Cryotubes, PP, sterile, Self-standing, 1.2 ml, Ext. diam. 12.6 mm, Height 40.2 mm, Thread Internal	100 pz
LLG09401160	Cryotubes, PP, sterile, Self-standing, 1.2 ml, Ext. diam. 12.6 mm, Height 40.2 mm, Thread External	100 pz.
LLG09401162	Cryotubes, PP, sterile, Self-standing, 2 ml, Ext. diam. 12.6 mm, Height 45.6 mm, Thread External	100 pz.
LLG09401166	Cryotubes, PP, sterile, Self-standing, 5 ml, Ext. diam. 12.6 mm, Height 87.6 mm, Thread External	100 pz
LLG09401164	Cryotubes, PP, sterile, Without ring, 2 ml, Ext. diam. 12.6 mm, Height 44.5 mm, Thread External	2 x 50 pz
LLG09401163	Cryotubes, PP, sterile, Self-standing, 2 ml, Ext. diam. 12.6 mm, Height 48.2 mm, Thread Internal	100 pezzi
LLG09401167	Cryotubes, PP, sterile, Without ring, 5 ml, Ext. diam. 12.6 mm, Height 89.2 mm, Thread Internal	100 pz.
LLG09401168	Cryotubes, PP, sterile, Self-standing, 5 ml, Ext. diam. 12.6 mm, Height 90.3 mm, Thread Internal	100 pz.
LLG09401165	Cryotubes, PP, sterile, Without ring, 2 ml, Ext. diam. 12.6 mm, Height 47 mm, Thread Internal	100 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

