



Racks for Cryogenic freezers CRYO 170 / CRYO 230

Racks for Cryogenic freezers CRYO 170 / CRYO 230



Code	Description	Packaging
LLG06267849	Racks for Cryogenic freezers CRYO 170 / CRYO 230, Type CFR100-06, dim. (W x D x H) 140 x 140 x 635 mm, Height Box 100 mm	1 pz.
LLG06281872	Racks for Cryogenic freezers CRYO 170 / CRYO 230, Type CFR75-05, dim. (W x D x H) 138 x 138 x 399 mm, Height Box 75 mm	
LLG06267848	Racks for Cryogenic freezers CRYO 170 / CRYO 230, Type CFR75-08, dim. (W x D x H) 140 x 140 x 629 mm, Height Box 75 mm	1 pz.
LLG06267847	Racks for Cryogenic freezers CRYO 170 / CRYO 230, Type CFR50-12, dim. (W x D x H) 140 x 140 x 665 mm, Height Box 50 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

