

**Orbital shaker LLG-uniSHAKER 25****LLG Labware****Orbital shaker LLG-uniSHAKER 25**

**Powerful and smooth running orbital shaker.** The maintenance-free brushless motor and the shaking mechanism with counterweight ensure smooth and stable running at all times. Digital display for speed and time. Two different platforms (610 x 460 mm or 760 x 460 mm) and 7 clamps for vessels from 25 ml to 2000 ml available. **Please order platforms and clamps separately.**

**-Maximum load 25 kg**

- Sturdy and robust steel housing
- Digital display with speed and time settings
- Variable speed from 50 to 500 min-1
- Orbital Diameter of 25 mm for thorough mixing
- Timer setting from 1 to 99 min and continuous mode
- Various platforms and brackets available

**Scope of supply:** shaker, mains cable

**Specifications**

**Type of movement:** orbital

**Orbital diameter:** 25 mm

**Speed range:**

**Max. load:** 50 ... 500 min-1

**25 kg**

**Motor:** brushless DC motor

**Timer:** 1 ... 99 min, continuous

**Power consumption:** 250 W

**Permissible ambient temperature:** 4 ... 75 °C

**Dimensions (W x D x H):** 600 x 550 x 150 mm

**Weight:** 48 kg

**Power supply:** 230 V, 50 Hz

**Warranty:** 3 years



| Code        | Description                                                                                        | Packaging |
|-------------|----------------------------------------------------------------------------------------------------|-----------|
| LLG06263589 | Accessories for Orbital shaker uniSHAKER25, Universal attachment, 610 x 460 mm, with securing bars | 1 pz.     |
| LLG06263591 | Accessories for Orbital shaker uniSHAKER25, Universal attachment, 610 x 460 mm                     | 1 pz.     |
| LLG06263592 | Accessories for Orbital shaker uniSHAKER25, Universal attachment, 760 x 460 mm                     | 1 pz.     |
| LLG06263593 | Accessories for Orbital shaker uniSHAKER25, Clamp, 100 ml                                          | 1 pz.     |
| LLG06263598 | Accessories for Orbital shaker uniSHAKER25, Adhesive rubber mat, blue, 150 x 150 x 3 mm            | 2 pz.     |
| LLG06263590 | uniSHAKER 25 orbital shaker, orbit 25 mm, 50-500 rpm, digital 230V with EU plug                    | 1 pz.     |
| LLG06263599 | uniSHAKER 25 Orbital shaker, orbit 25 mm, 50-500 rpm, digital 230V, with UK plug                   | 1 pz.     |
| LLG06263594 | Accessories for Orbital shaker uniSHAKER25, Clamp, 250 ml                                          | 1 pz.     |
| LLG06263595 | Accessories for Orbital shaker uniSHAKER25, Clamp, 500 ml                                          | 1 pz.     |
| LLG06263596 | Accessories for Orbital shaker uniSHAKER25, Clamp, 1000 ml                                         | 1 pz.     |
| LLG06263597 | Accessories for Orbital shaker uniSHAKER25, Clamp, 2000 ml                                         | 1 pz.     |
| LLG06263699 | Accessories for Orbital shaker uniSHAKER25, Clamp, 25 ml                                           | 1 pz.     |
| LLG06263700 | Accessories for Orbital shaker uniSHAKER25, Clamp, 50 ml                                           | 1 pz.     |





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

