

**Orbital shaker LLG-uniSHAKER 25****LLG Labware****Orbital shaker LLG-uniSHAKER 25**

Powerful and smooth running orbital shaker. The maintenance-free brushless motor and the shaking mechanism with counterweight ensure smooth and stable running at all times. Digital display for speed and time. Two different platforms (610 x 460 mm or 760 x 460 mm) and 7 clamps for vessels from 25 ml to 2000 ml available. **Please order platforms and clamps separately.**

-Maximum load 25 kg

- Sturdy and robust steel housing
- Digital display with speed and time settings
- Variable speed from 50 to 500 min⁻¹
- Orbital Diameter of 25 mm for thorough mixing
- Timer setting from 1 to 99 min and continuous mode
- Various platforms and brackets available

Scope of supply: shaker, mains cable

Specifications

Type of movement: orbital

Orbital diameter: 25 mm

Speed range:

Max. load: 50 ... 500 min⁻¹
25 kg

Motor: brushless DC motor

Timer: 1 ... 99 min, continuous

Power consumption: 250 W

Permissible ambient temperature: 4 ... 75 °C

Dimensions (W x D x H): 600 x 550 x 150 mm

Weight: 48 kg

Power supply: 230 V, 50 Hz

Warranty: 3 years



Code	Description	Packaging
LLG06263589	Accessories for Orbital shaker uniSHAKER25, Universal attachment, 610 x 460 mm, with securing bars	1 pz.
LLG06263591	Accessories for Orbital shaker uniSHAKER25, Universal attachment, 610 x 460 mm	1 pz.
LLG06263592	Accessories for Orbital shaker uniSHAKER25, Universal attachment, 760 x 460 mm	1 pz.
LLG06263593	Accessories for Orbital shaker uniSHAKER25, Clamp, 100 ml	1 pz.
LLG06263598	Accessories for Orbital shaker uniSHAKER25, Adhesive rubber mat, blue, 150 x 150 x 3 mm	2 pz.
LLG06263590	uniSHAKER 25 orbital shaker, orbit 25 mm, 50-500 rpm, digital 230V with EU plug	1 pz.
LLG06263599	uniSHAKER 25 Orbital shaker, orbit 25 mm, 50-500 rpm, digital 230V, with UK plug	1 pz.
LLG06263594	Accessories for Orbital shaker uniSHAKER25, Clamp, 250 ml	1 pz.
LLG06263595	Accessories for Orbital shaker uniSHAKER25, Clamp, 500 ml	1 pz.
LLG06263596	Accessories for Orbital shaker uniSHAKER25, Clamp, 1000 ml	1 pz.
LLG06263597	Accessories for Orbital shaker uniSHAKER25, Clamp, 2000 ml	1 pz.
LLG06263699	Accessories for Orbital shaker uniSHAKER25, Clamp, 25 ml	1 pz.
LLG06263700	Accessories for Orbital shaker uniSHAKER25, Clamp, 50 ml	1 pz.





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



*CARLO ERBA
Reagents operates with
a Certified Quality
Management System*

