

**Sampling Set with Spatulas SteriPlast Kit****Bürkle****Sampling Set with Spatulas SteriPlast® Kit**

The SteriPlast® Kit is the ideal combination for sampling where purity and sterility is essential. All parts are gamma-sterilized and suitable for use with foodstuffs. The spatulas are Cleanroom manufactured and individually packaged. The sample bag SteriBag Premium has an integrated sealing and is liquid-tight. An unique tracking code on both the bag and the tear-off field guarantees complete documentation and traceability.

Scope of supply: 10 x SteriPlast® spatulas 150 mm insertion depth, and 10 x SteriBag Premium bags 300 ml.



Code	Description	Packaging
LLG06283733	Sampling Set with Spatulas SteriPlast® Kit, Type SteriPlast® Kit (Spatulas and Bags)	1 kit



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

