

**Sampling Set with Spatulas SteriPlast Kit****Bürkle****Sampling Set with Spatulas SteriPlast® Kit**

The SteriPlast® Kit is the ideal combination for sampling where purity and sterility is essential. All parts are gamma-sterilized and suitable for use with foodstuffs. The spatulas are Cleanroom manufactured and individually packaged. The sample bag SteriBag Premium has an integrated sealing and is liquid-tight. An unique tracking code on both the bag and the tear-off field guarantees complete documentation and traceability.

**Scope of supply:** 10 x SteriPlast® spatulas 150 mm insertion depth, and 10 x SteriBag Premium bags 300 ml.



Code	Description	Packaging
LLG06283733	Sampling Set with Spatulas SteriPlast® Kit, Type SteriPlast® Kit (Spatulas and Bags)	1 kit



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

