

**Sampling Set SteriPlast Kit, sterile****Bürkle****Sampling Set SteriPlast Kit, sterile**

The SteriPlast® Kit is the ideal combination for sampling where purity and sterility is essential. All parts are sterilised by gamma rays and suitable for use with foodstuffs. The scoops or spatulas are Cleanroom manufactured and individually packaged. The tamper-proof sample bag SteriBag Premium has an integrated sealing and is liquid-tight. A unique tracking code on both the bag and the tear-off field guarantees complete documentation and traceability.

Scope of supply: 10 x SteriPlast scoops 50 ml and 10 x SteriBag Premium bags 300 ml.



Code	Description	Packaging
LLG06283732	Sampling Set SteriPlast Kit, sterile, Type SteriPlast Kit (Scoops and Bags)	1 kit



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

