



## Analysis- and precision balances Secura

### Analysis- and precision balances Secura®

Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, Secura® also features built-in protection systems for complete reliability and regulatory compliance, the safe and secure way. Real-time guidance prompts for leveling, automatic internal adjustment and 100 % traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Intuitive Operation
- Ergonomic Draft Shield
- Automatic Internal Adjustment
- Real-Time Level Support
- Type CEU is factory veriflicated (EU)

**Specifications Secura® 124/213/224/313/513/613 |  
612/1102/2102/3102/5102/6102**

**Overall dimensions (W x D x H) 360 x 216 x 320 mm | 360 x 216 x 95 mm**

**Weight: 5.1 kg | 5.2 kg**

**Height weighing space: 209 mm | -**

**Power supply: 100-240 V, 50/60 Hz**



Code	Description	Packaging
LLG06262434	Analysis- and precision balances Secura®, Type SECURA124-1S, 120 g	1 pz.
LLG06282260	Analysis- and precision balances Secura®, Type SECURA125-1S, 120 g	1 pz.
LLG06262432	Analysis- and precision balances Secura®, Type SECURA224-1S, 220 g	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

