



Analysis- and precision balances Secura

Analysis- and precision balances Secura®

Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, Secura® also features built-in protection systems for complete reliability and regulatory compliance, the safe and secure way. Real-time guidance prompts for leveling, automatic internal adjustment and 100 % traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Intuitive Operation
- Ergonomic Draft Shield
- Automatic Internal Adjustment
- Real-Time Level Support
- Type CEU is factory veriflicated (EU)


Specifications Secura® 124/213/224/313/513/613 | 612/1102/2102/3102/5102/6102

Overall dimensions (W x D x H) 360 x 216 x 320 mm | 360 x 216 x 95 mm

Weight: 5.1 kg | 5.2 kg

Height weighing space: 209 mm | -

Power supply: 100-240 V, 50/60 Hz



Code	Description	Packaging
LLG06262434	Analysis- and precision balances Secura®, Type SECURA124-1S, 120 g	1 pz.
LLG06282260	Analysis- and precision balances Secura®, Type SECURA125-1S, 120 g	1 pz.
LLG06262432	Analysis- and precision balances Secura®, Type SECURA224-1S, 220 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

