



DURAN Filtering Flask with Glass Hose Connection, Erlenmeyer shape

DWK Life Sciences GmbH
(Duran)

DURAN® Filtering Flask with Glass Hose Connection, Erlenmeyer shape

Due to the heavy wall thickness the apparatus is vacuum-tight.



Code	Description	Packaging
LLG09051201	Filtering flasks 100ml with glass olive	1 pz.
LLG09051202	Filtering flasks 250ml with glass olive	1 pz.
LLG09051210	Filtering flasks 1L with glass olive	1 pz
LLG09051220	Filtering flasks 2L with glass olive	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

