

**Digital Vortexer LLG-uniTEXER 4****LLG Labware****Digital Vortexer LLG-uniTEXER 4**

Compact vortexer with brushless motor and infinitely variable speed control up to 4200 rpm for optimal mixing results. The digital display alternately shows the set speed and time.

- Programmable Pulse mode
- 3-way switch (on/off/touch)
- Timer adjustable from 1 min up to 999 min

Scope of supply: Vortexer, conical tube attachment, round disk attachment, universal plug

Specifications :

Type of movement: Orbital

Orbital diameter: 4 mm

Speed range: 300 ... 4200 rpm

Dimensions (W x D x H): 205 x 136 x 138.5 mm

Weight: 2.8 kg

Power supply: 110 ... 240 V, 50/60 Hz

IP code: IP 42

Warranty: 3 years



Code	Description	Packaging
LLG06263500	Digital Vortexer LLG-uniTEXER 4	1 pc



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

