

**Rotavapor R-300 with heating bath B-301****BÜCHI****Rotavapor® R-300 with heating bath B-301**

The rotary evaporator R-300 meets the highest needs in convenience, productivity and versatility. The rotation speed, heating temperature and lift position is shown digitally on the integrated display. Automatic change of the direction of rotation programmable.

**Lift System**

Option 1: Manual lift

Option 2: Electronic lift

**Heating Bath**

**Option 1: B-301.** Temperature range of 20 °C to 95 °C (water). Max. evaporating flask size: 1000 mL (3 kg). Cordless connection. Digital display of set and actual temperature, rotation speed, lift position. Lock of heating temperature. Fast heat up time.

**Option 2: B-305.** Temperature range of 20 °C to 220 °C. Max. evaporating flask size: 5000 mL (3 kg). Cordless connection. Digital display of set and actual temperature, rotation speed, lift position. Lock of heating temperature.

**Standard Joint**

SJ 29/32

**Glass Assembly**

**Option V: Vertical condenser (V):** Vertical condenser used along with a recirculating chiller or tap water.

**Option C: Cold trap (C):** Cold trap used for low boiling solvents along with e.g. dry ice for cooling.

Other glass assembly configuration on request

**Protective Coating**

Option: Transparent PLASTIC+GLASS (P+G) safety coating. Protect against mechanical damage. In the event of breakage there is no loss of sample.



Code	Description	Packaging
LLG06267887	Rotavapor® R-300 with heating bath B-301, Type R-300, Lift manually, Capacitor vertical, Coating no	1 pz.
LLG06267889	Rotavapor® R-300 with heating bath B-301, Type R-300, Lift electric, Capacitor cold trap, Coating yes	1 pz.
LLG06267888	Rotavapor® R-300 with heating bath B-301, Type R-300, Lift electric, Capacitor vertical, Coating yes	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

