

**Beakers glass, DURAN, low form****DWK Life Sciences****Beakers glass, DURAN®, low form**

DIN 12331. ISO 3819. With approximate graduations and spout. With easy-to-read scale and large labelling field for easy marking. In fired-on, highly durable white ceramic. Due to the uniform wall thickness, it is suitable for use as a heating vessel. Glass Type I / neutral as per USP, EP and JP. With retrace code (Batch Identification, for all) except 5ml, 10ml and 25ml. Certificate available via the Internet. Autoclavable.



Code	Description	Packaging
LLG09013053	Beakers glass, DURAN®, low form, Nominal capacity 800 ml	1 pz
LLG09013054	Beakers glass, DURAN®, low form, Nominal capacity 1000 ml	1 pz
LLG09013063	Beakers glass, DURAN®, low form, Nominal capacity 2000 ml	1 pz
LLG09013068	Beakers glass, DURAN®, low form, Nominal capacity 3000 ml	1 pz
LLG09013073	Beakers glass, DURAN®, low form, Nominal capacity 5000 ml	1 pz
LLG09013086	Beakers glass, DURAN®, low form, Nominal capacity 10000 ml	1 pz
LLG09013007	Beakers glass, DURAN®, low form, Nominal capacity 5 ml	1 pz
LLG09013008	Beakers glass, DURAN®, low form, Nominal capacity 10 ml	1 pz
LLG09013014	Beakers glass, DURAN®, low form, Nominal capacity 25 ml	1 pz
LLG09013017	Beakers glass, DURAN®, low form, Nominal capacity 50 ml	1 pz
LLG09013024	Beakers glass, DURAN®, low form, Nominal capacity 100 ml	1 pz
LLG09013029	Beakers glass, DURAN®, low form, Nominal capacity 150 ml	1 pz
LLG09013036	Beakers glass, DURAN®, low form, Nominal capacity 250 ml	1 pz
LLG09013041	Beakers glass, DURAN®, low form, Nominal capacity 400 ml	1 pz
LLG09013048	Beakers glass, DURAN®, low form, Nominal capacity 600 ml	1 pz





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

