

**LLG-Cell strainers, Nylon, sterile****LLG Labware****LLG-Cell strainers, Nylon, sterile**

Cell strainers are especially suited for the generation of homogenous single cell suspensions from tissue or primary cell lysates

- Available in 3 mesh sizes, 40 µm, 70 µm, and 100 µm
- Colour coded mesh size
- Designed to fit into 50 ml tubes
- Individually packaged
- Handle and mantle surface for improved aseptic handling
- Nylon mesh
- Pyrogen-free



Code	Description	Packaging
LLG06268093	Cell strainers, Nylon, sterile, Mesh size 40 µm, Colour code blue	50 pezzi
LLG06268092	Cell strainers, Nylon, sterile, Mesh size 70 µm, Colour code white	50 pz
LLG06268091	Cell strainers, Nylon, sterile, Mesh size 100 µm, Colour code yellow	50 pezzi



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

