



Accessories for Rotor Beater Mill SR 300

Accessories for Rotor Beater Mill SR 300

- Ring sieves 360° with trapezoid holes, stainless
- Retaining frame 360°, stainless

Further aperture sizes on request.



Code	Description	Packaging
LLG09738914	Accessories for Rotor Beater Mill SR 300, Ring sieve, Aperture size 0,25 mm	1 pz.
LLG09738915	Accessories for Rotor Beater Mill SR 300, Ring sieve, Aperture size 0,50 mm	1 pz.
LLG09738916	Accessories for Rotor Beater Mill SR 300, Ring sieve, Aperture size 0,75 mm	1 pz.
LLG09738917	Accessories for Rotor Beater Mill SR 300, Ring sieve, Aperture size 1,00 mm	1 pz.
LLG06268107	Retaining frame	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

