

**Beakers glass, DURAN, tall form****DWK Life Sciences****Beakers glass, DURAN®, tall form**

DIN 12331. ISO 3819. With approximate graduations and spout. With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Due to the uniformly wall thickness, it is suitable for use as a heating vessel. Glass Type I / neutral glass as per USP, EP and JP. With retrace code (Batch Identification). Certificate available via the Internet. Autoclavable.



| Code        | Description  | Packaging |
|-------------|--|-----------|
| LLG09013163 | Beakers glass, DURAN®, tall form, Nominal capacity 2000 ml | 1 pz      |
| LLG09013168 | Beakers glass, DURAN®, tall form, Nominal capacity 3000 ml | 1 pz      |
| LLG09013117 | Beakers glass, DURAN®, tall form, Nominal capacity 50 ml   | 1 pz      |
| LLG09013124 | Beakers glass, DURAN®, tall form, Nominal capacity 100 ml  | 1 pz      |
| LLG09013129 | Beakers glass, DURAN®, tall form, Nominal capacity 150 ml  | 1 pz      |
| LLG09013136 | Beakers glass, DURAN®, tall form, Nominal capacity 250 ml  | 1 pz      |
| LLG09013141 | Beakers glass, DURAN®, tall form, Nominal capacity 400 ml  | 1 pz      |
| LLG09013148 | Beakers glass, DURAN®, tall form, Nominal capacity 600 ml  | 1 pz      |
| LLG09013153 | Beakers glass, DURAN®, tall form, Nominal capacity 800 ml  | 1 pz      |
| LLG09013154 | Beakers glass, DURAN®, tall form, Nominal capacity 1000 ml | 1 pz      |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

