

**Beakers glass, DURAN, tall form****DWK Life Sciences****Beakers glass, DURAN®, tall form**

DIN 12331. ISO 3819. With approximate graduations and spout. With easy-to-read scale and large labelling field for easy marking, in fired-on, highly durable white ceramic. Due to the uniformly wall thickness, it is suitable for use as a heating vessel. Glass Type I / neutral glass as per USP, EP and JP. With retrace code (Batch Identification). Certificate available via the Internet. Autoclavable.



Code	Description	Packaging
LLG09013163	Beakers glass, DURAN®, tall form, Nominal capacity 2000 ml	1 pz
LLG09013168	Beakers glass, DURAN®, tall form, Nominal capacity 3000 ml	1 pz
LLG09013117	Beakers glass, DURAN®, tall form, Nominal capacity 50 ml	1 pz
LLG09013124	Beakers glass, DURAN®, tall form, Nominal capacity 100 ml	1 pz
LLG09013129	Beakers glass, DURAN®, tall form, Nominal capacity 150 ml	1 pz
LLG09013136	Beakers glass, DURAN®, tall form, Nominal capacity 250 ml	1 pz
LLG09013141	Beakers glass, DURAN®, tall form, Nominal capacity 400 ml	1 pz
LLG09013148	Beakers glass, DURAN®, tall form, Nominal capacity 600 ml	1 pz
LLG09013153	Beakers glass, DURAN®, tall form, Nominal capacity 800 ml	1 pz
LLG09013154	Beakers glass, DURAN®, tall form, Nominal capacity 1000 ml	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

