



Accessories for Overhead Stirrers Hei-TORQUE

Accessories for Overhead Stirrers Hei-TORQUE



Code	Description	Packaging
LLG06228500	Accessories for Overhead Stirrers Hei-TORQUE, Type Base Stand S2 XXL, For Hei-TORQUE	1 pz.
LLG07047043	Accessories for Overhead Stirrers Hei-TORQUE, Type Flexible shaft including chuck, 8mm, For Hei-TORQUE	1 pz.
LLG07621811	Accessories for Overhead Stirrers Hei-TORQUE, Type Telescopic stand, For Hei-TORQUE	1 pz.
LLG09816566	Accessories for Overhead Stirrers Hei-TORQUE, Type Stirrer shaft guard, PMMA, height-adjustable, For Hei-TORQUE Core	1 pz.
LLG09816508	Chuck 8,5mm span	1 pz.
LLG09816563	Chuck, clamping range 10.5 mm	1 pz.
LLG09816562	Accessories for Overhead Stirrers Hei-TORQUE, Type Stand clamp for holder rod 13 to 32mm, For Hei-TORQUE	1 pz.
LLG09816564	Accessories for Overhead Stirrers Hei-TORQUE, Type Stand S 2, 425 x 420 x 700mm, 5.8kg, For Hei-TORQUE	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

