



## Accessories for Overhead Stirrers Hei-TORQUE

### Accessories for Overhead Stirrers Hei-TORQUE



| Code        | Description                                                                                                          | Packaging |
|-------------|----------------------------------------------------------------------------------------------------------------------|-----------|
| LLG06228500 | Accessories for Overhead Stirrers Hei-TORQUE, Type Base Stand S2 XXL, For Hei-TORQUE                                 | 1 pz.     |
| LLG07047043 | Accessories for Overhead Stirrers Hei-TORQUE, Type Flexible shaft including chuck, 8mm, For Hei-TORQUE               | 1 pz.     |
| LLG07621811 | Accessories for Overhead Stirrers Hei-TORQUE, Type Telescopic stand, For Hei-TORQUE                                  | 1 pz.     |
| LLG09816566 | Accessories for Overhead Stirrers Hei-TORQUE, Type Stirrer shaft guard, PMMA, height-adjustable, For Hei-TORQUE Core | 1 pz.     |
| LLG09816508 | Chuck 8,5mm span                                                                                                     | 1 pz.     |
| LLG09816563 | Chuck, clamping range 10.5 mm                                                                                        | 1 pz.     |
| LLG09816562 | Accessories for Overhead Stirrers Hei-TORQUE, Type Stand clamp for holder rod 13 to 32mm, For Hei-TORQUE             | 1 pz.     |
| LLG09816564 | Accessories for Overhead Stirrers Hei-TORQUE, Type Stand S 2, 425 x 420 x 700mm, 5.8kg, For Hei-TORQUE               | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

