



Deionisers, pressure-resistant, behropur

Deionisers, pressure-resistant, behropur®

Compact and secure mixed bed deionisers for small to medium-sized volumes of ultra pure water. Ideal for feeding laboratory washing machines, general requirements in the laboratory, low level consumption in industry and for post-treatment desalination of the output from reverse osmosis systems.

- Optimal utilisation of the deioniser due to totally uniform water distribution in the resin chamber.
- Can be connected directly to the raw water mains without a pressure reducer.
- Back pressure resistant.
- Also available with water quality cut-out and solenoid valve, which engages if specified limits are exceeded, and reservoir level control.
- Resistivity control available directly on the unit or remotely at any location.

Performance data: B10dN / B22dN / B45dN

Cation exchange capacity at 10°dH:

Flow max.:

Dia.:

Height incl. LF: 1200 L / 2400 L / 5500 L

300 / 500 / 800 L/hr

21 / 21 / 26 cm

68 / 112 / 125 cm

Height, cartridge only: 55 / 98 / 110 cm



Code	Description	Packaging
LLG09882123	Deionisers, pressure-resistant, behropur®, Type B10dNZ, behropur water de-ioniser cartridge, nylon, pressure resistant for 9.882 120	1 pz.
LLG09882124	Deionisers, pressure-resistant, behropur®, Type B22dNZ, behropur water de-ioniser cartridge, nylon, pressure resistant for 9.882 121	1 pz.
LLG09882125	Deionisers, pressure-resistant, behropur®, Type B45dNZ, behropur water de-ioniser cartridge, nylon, pressure resistant for 9.882 122	1 pz.
LLG09882126	Deionisers, pressure-resistant, behropur®, Type B10dNA, Pressure-resistant mixed bed unit, with resistivity meter	1 pz.
LLG09882128	Deionisers, pressure-resistant, behropur®, Type B45dNA, Pressure-resistant mixed bed unit, meter with water quality cut-out	1 pz.
LLG09882127	Deionisers, pressure-resistant, behropur®, Type B22dNA, Pressure-resistant mixed bed unit stivity, meter with water quality cut-out	1 pz.
LLG09882120	Deionisers, pressure-resistant, behropur®, Type B10dN, Pressure-resistant mixed bed unit, with resistivity meter	1 pz.
LLG09882121	Deionisers, pressure-resistant, behropur®, Type B22dN, Pressure-resistant mixed bed unit, with resistivity meter	1 pz.
LLG09882122	Deionisers, pressure-resistant, behropur®, Type B45dN, Pressure-resistant mixed bed unit, with resistivity meter	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



