

**LLG-Transfer membrane, nitrocellulose****LLG Labware****LLG-Transfer membrane, nitrocellulose**

Pure Nitrocellulose transfer membrane for all protein or immunoblotting applications. The high sensitivity and low background ensures excellent results in all transfers. Compatible with all detection systems. BSA binding capacity up to 100 µg/cm<sup>2</sup>.



Code	Description	Packaging
LLG06264061	Transfer membrane, nitrocellulose, Width 300 mm, Length 3000 mm, Pore size 0.22 µm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

