

**Syringeless filter system Claristep**

Sartorius Lab Instruments

**Syringeless filter system Claristep®**

Sartorius has developed a easy-to-use and straightforward filtration setup. The manually operated Claristep® Filtration System offers a way for clarifying your samples prior to analysis.

- Up to 8 samples are processed simultaneously
- No syringe required
- No need for a vacuum source or a power supply
- For low sample volumes ranging from 60 to 600 µl
- Hold-up volume 30 µl
- Suitable for all 12 x 32 mm test tubes

**Scope of supply:** Claristep® Station consisting of a base, a lid and an exchangeable tray



Code	Description	Packaging
LLG06283802	Syringeless filter system Claristep®, Type Claristep® Single Tray (replacement)	1 pz.
LLG06283801	Syringeless filter system Claristep®, Type Claristep® Station complete	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

