

**Filters for Syringeless filter system Claristep**

Sartorius Lab Instruments

Filters for Syringeless filter system Claristep®

Housing made of polypropylene (PP), membrane made of regenerated cellulose (RC). Effective filtration diameter 9.7 mm.



Code	Description	Packaging
LLG06283797	Filters for Syringeless filter system Claristep®, Pore size 0.2 µm, Membrane RC	96 pz.
LLG06283798	Filters for Syringeless filter system Claristep®, Pore size 0.2 µm, Membrane RC	480 pz.
LLG06283799	Filters for Syringeless filter system Claristep®, Pore size 0.45 µm, Membrane RC	96 pz.
LLG06283800	Filters for Syringeless filter system Claristep®, Pore size 0.45 µm, Membrane RC	480 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

