

**Griffin beakers, ETFE****VITLAB****Griffin beakers, ETFE**

Transparent. With easily readable, printed black scale.

Very good chemical resistance and very high thermal stability from -100°C to +150°C.

To preserve markings, cleaning at greater than 60°C is not recommended.



Code	Description	Packaging
LLG09013704	Griffin beakers, ETFE, 100 ml, Grad. 20 ml, Diam. 50 mm, Height 72 mm	1 pz.
LLG09013710	Griffin beakers, ETFE, 1000 ml, Grad. 100 ml, Diam. 105 mm, Height 143 mm	1 pz.
LLG09013708	Griffin beakers, ETFE, 600 ml, Grad. 100 ml, Diam. 91 mm, Height 125 mm	1 pz.
LLG09013702	Griffin beakers, ETFE, 25 ml, Grad. 5 ml, Diam. 32 mm, Height 50 mm	1 pz.
LLG09013706	Griffin beakers, ETFE, 250 ml, Grad. 50 ml, Diam. 67 mm, Height 96 mm	1 pz.
LLG09013707	Griffin beakers, ETFE, 400 ml, Grad. 50 ml, Diam. 77 mm, Height 109 mm	1 pz.
LLG09013703	Griffin beakers, ETFE, 50 ml, Grad. 10 ml, Diam. 39 mm, Height 59 mm	1 pz.
LLG09013192	Griffin beakers, ETFE, 500 ml, Grad. 100 ml, Diam. 88 mm, Height 122 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

