



Laboratory Max./Min. thermometer

Laboratory Max./Min. thermometer

Display of current, highest and lowest temperature, hold function, alarm at freely adjustable temperature values optical and audibly (upper/lower limit), selectable °C/°F readout. Supplied with table stand and magnet. 3V Lithium CR2032 button cell (acc. to EN 13485).

Specifications LT101 | LT102

Sensor:

Resolution: NTC

0.1 °C

Response time t90: appr. 15 seconds in water

Dimensions (L x W x H): 86 x 57 x 30 mm

Sensor tube (L x diam.): 100 x 3.5 mm | 30 x 5 mm

Weight: 100 g



Code	Description	Packaging
LLG09851022	Laboratory Max./Min. thermometer, Type LT102	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

