



## LLG Syringe Filters SPHEROS, PTFE

### LLG Syringe Filters *SPHEROS*, PTFE

Cost-effective syringe filters for filtration of a broad variety of solvents and aqueous or inorganic solutions. These syringe filters cover most applications in HPLC, pharmaceutical, environmental, biotechnology, and food and beverage testing laboratories.

- Robust Polypropylene housing
  - Multifunctional syringe filters with female Luer-Lock inlet and male Luer-Slip outlet
  - Suitable for all syringes with Luer connection
  - Sterile or non-sterile
  - Sterile products in individual hard blister packaging and dispenser box
- Hydrophobic membrane for filtration of non-polar liquids and gases. Very good chemical resistance to all kinds of solvents and alkaline solutions. With female Luer-Lock inlet and male Luer-Slip outlet.



Code	Description	Packaging
LLG06272816	LLG Syringe FiltersSPHEROS, PTFE, Filter diam. 13 mm, Pore size 0.22 µm, Housing PP, white, Sterile -	100 pz.
LLG06272817	LLG Syringe FiltersSPHEROS, PTFE, Filter diam. 13 mm, Pore size 0.45 µm, Housing PP, white, Sterile -	100 pz.
LLG06272818	LLG Syringe FiltersSPHEROS, PTFE, Filter diam. 25 mm, Pore size 0.22 µm, Housing PP, white, Sterile -	50 pz.
LLG06272819	LLG Syringe FiltersSPHEROS, PTFE, Filter diam. 25 mm, Pore size 0.45 µm, Housing PP, white, Sterile -	50 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

