

**Analytical balances Adventurer AX****OHAUS****Analytical balances Adventurer® AX**

- Fast stabilization and reliable operation to ensure solid weighing performance and accurate, repeatable results
- Two USB ports with GLP/GMP tracking capabilities assist the monitoring and reporting of the data
- The wide viewing angle colour touchscreen provides access to all applications
- Reduced footprint due to two pieces, top-mounted draft shield doors
- AutoCalä automatic internal calibration system or traditional external calibration (external calibration models -/E)
- EC Type Approved models available
- Applications: Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, dynamic weighing, totalization/statistics, formulation, density determination, display hold
- Display: 4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness
- Operation: Mains operation
- Communication: RS232, USB Host, USB Device (included). GLP/GMP compliant data output
- Construction: Die-cast metal base, ABS top housing, stainless steel pan, glass draft shield with two piece top mounted side doors and sliding top door, replaceable in-use cover
- Design features: Two-piece, top mounted, space-saving draftshield with removable panels on 0.1mg and 0.001g models. AutoCalä automatic internal calibration or external calibration. Level assist feature with illuminated level indicator. Selectable environmental and auto-print settings. Integral weigh below hook. Protective full housing in-use cover. Software lock switch

**Specifications****Readability: 0.1 mg****Ø weighing plate: 90 mm****Power supply: 100...240 V/50/60 Hz**

| Code | Description | Packaging |
|-------------|--|-----------|
| LLG06268193 | Analytical balances Adventurer® AX, Type AX124, Weighing range 120 g | 1 pz. |
| LLG06268195 | Analytical balances Adventurer® AX, Type AX124/E, Weighing range 120 g | 1 pz. |
| LLG06268194 | Analytical balances Adventurer® AX, Type AX124M, Weighing range 120 g | |
| LLG06268196 | Analytical balances Adventurer® AX, Type AX224, Weighing range 220 g | 1 pz. |
| LLG06268198 | Analytical balances Adventurer® AX, Type AX224/E, Weighing range 220 g | 1 pz. |
| LLG06268197 | Analytical balances Adventurer® AX, Type AX224M, Weighing range 220 g | 1 pz. |
| LLG06268199 | Analytical balances Adventurer® AX, Type AX324, Weighing range 320 g | 1 pz. |
| LLG06268200 | Analytical balances Adventurer® AX, Type AX324M, Weighing range 320 g | |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

