

**Precision Balances Adventurer AX****OHAUS****Precision Balances Adventurer® AX**

- Fast stabilization and reliable operation to ensure solid weighing performance and accurate, repeatable results
- Two USB ports with GLP/GMP tracking capabilities assist the monitoring and reporting of the data
- The wide viewing angle 4.3" colour touchscreen provides access to all applications
- Reduced footprint due to two piece, top-mounted draft shield doors
- The 0.01g models feature large weighing pans
- AutoCalä automatic internal calibration system or traditional external calibration (external calibration models -/E)
- EC Type Approved models available
- Applications: Weighing (17 units + custom unit), parts counting, percentage weighing, check weighing, dynamic weighing, totalization/statistics, formulation, density determination, display hold
- Display: 4.3" (109 mm) Full-colour VGA graphic touchscreen with user-controlled brightness
- Operation: Mains operation
- Communication: RS232, USB Host, USB Device (included). GLP/GMP compliant data output with real-time clock
- Construction: Die-cast metal base, ABS top housing, stainless steel pan, glass draft shield with two piece top mounted side doors and sliding top door (0.1 and 1mg models only), illuminated up-front level indicator, integral weigh below hook, security bracket, calibration lock, full housing in-use cover
- Design features: Two-piece, top mounted, space-saving draft shield with removable panels on 0.1mg and 0.001g models. AutoCalä automatic internal calibration or external calibration. Level assist feature with illuminated level indicator. Selectable environmental and auto-print settings. Integral weigh below hook. Protective full housing in-use cover. Software lock switch



Code	Description	Packaging
LLG06267263	Precision Balances Adventurer® AX, Type AX2202, Weighing range 2200 g	1 pz.
LLG06268219	Precision Balances Adventurer® AX, Type AX1502M, Weighing range 1520 g	
LLG06268220	Precision Balances Adventurer® AX, Type AX1502/E, Weighing range 1520 g	1 pz.
LLG06268222	Precision Balances Adventurer® AX, Type AX2202/E, Weighing range 2200 g	1 pz.
LLG06268224	Precision Balances Adventurer® AX, Type AX4202M, Weighing range 4200 g	
LLG06268225	Precision Balances Adventurer® AX, Type AX4202/E, Weighing range 4200 g	1 pz.
LLG06268226	Precision Balances Adventurer® AX, Type AX5202, Weighing range 5200 g	1 pz.
LLG06280937	Precision Balances Adventurer® AX, Type AX422, Weighing range 420 g	1 pz.
LLG06280940	Precision Balances Adventurer® AX, Type AX822/E, Weighing range 820 g	1 pz.
LLG06280941	Precision Balances Adventurer® AX, Type AX2201, Weighing range 2200 g	1 pz.
LLG06280942	Precision Balances Adventurer® AX, Type AX2201/E, Weighing range 2200 g	1 pz.
LLG06280943	Precision Balances Adventurer® AX, Type AX4201, Weighing range 4200 g	1 pz.
LLG06280944	Precision Balances Adventurer® AX, Type AX4201/E, Weighing range 4200 g	1 pz.
LLG06280946	Precision Balances Adventurer® AX, Type AX8201/E, Weighing range 8200 g	1 pz.
LLG06280947	Precision Balances Adventurer® AX, Type AX8201/M, Weighing range 8200 g	
LLG06268209	Precision Balances Adventurer® AX, Type AX223M, Weighing range 220 g	
LLG06268210	Precision Balances Adventurer® AX, Type AX423, Weighing range 420 g	1 pz.
LLG06268211	Precision Balances Adventurer® AX, Type AX423M, Weighing range 420 g	
LLG06268212	Precision Balances Adventurer® AX, Type AX423/E, Weighing range 420 g	1 pz.
LLG06268213	Precision Balances Adventurer® AX, Type AX523, Weighing range 520 g	1 pz.



Code	Description	Packaging
LLG06268214	Precision Balances Adventurer® AX, Type AX523M, Weighing range 520 g	
LLG06268215	Precision Balances Adventurer® AX, Type AX523/E, Weighing range 520 g	1 pz.
LLG06268218	Precision Balances Adventurer® AX, Type AX1502, Weighing range 1520 g	1 pz.
LLG06268208	Precision Balances Adventurer® AX, Type AX223, Weighing range 220 g	1 pz.
LLG06266893	Precision Balances Adventurer® AX, Type AX223/E, Weighing range 220 g	1 pz.
LLG06268223	Precision Balances Adventurer® AX, Type AX4202, Weighing range 4200 g	1 pz.
LLG06268227	Precision Balances Adventurer® AX, Type AX5202M, Weighing range 5200 g	
LLG06280945	Precision Balances Adventurer® AX, Type AX8201, Weighing range 8200 g	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

