

**Universal shakers, SM 30 series / SM 30 control****Bühler****Universal shakers, SM 30 series / SM 30 control**

Powerful shakers for large loads. Choice of orbital or reciprocating action models. Suitable for continuous operation. For analytical or diagnostic processes in biology, chemistry or pharmacy, plant culture or in the food industry, e.g. for transport simulations. Various attachments for different requirements are available as accessories. Programmable controller with Touch-Display and advanced features for control devices.

**Scope of supply:** supplied with shaking rack and rubber mat.

Attachments required, please order separately.

Speed range: 15 ... 300 rpm

stepless, electronic adjustment, control devices in steps of 5

Orbit or travel: 26 mm (A + C models)

30, 46 or 50 mm (B models)

Max. loading: 30 kg

Timer: 0 ... 120 min / continuous, control devices programmable / continuous

Platform (W x D): 560 x 400 mm

Overall (W x D x H): 680 x 610 x 160 mm

Weight: 33 kg

Supply requirements: 230 V, 50/60 Hz



Code	Description	Packaging
LLG09838101	Universal shakers SM 30A,orbital moti/1	1 pz.
LLG09838115	Universal shakers SM 30C,orbital/reci/1	1 pz.
LLG09838135	Universal shakers SM 30C-control,with/1	1 pz.
LLG09838111	Universal shakers SM 30,B,reciprocat./1	1 pz.
LLG09838131	Universal shakers SM 30B-control,reci/1	1 pz.
LLG09838121	Universal shakers SM 30A-control,orbi/1	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

