

**Accessories for universal shaker SM 30****Bühler****Accessories for universal shaker SM 30**

Combifix SM A attachment: Flexible attachment with 5 clamping bars for various containers with flat bases.

Combifix SM B attachment: With 4 clamping bars for holding horizontal or tall containers securely.

Combifix SM C attachment: For separating funnels.

Universal tray SM: For flexible fitting with accessory spring clips or test tube stands.

Sliding Plates: Sliding plates with Combifix SM (rack systems Combifix SM, but with additional base plate).

Mounting set: for universal trays, or sliding plates on the 2-storey top frame SM.

SM attachment, 2-storey support: With 2 boards for spring clips (please order separately) or for mounting Combifix SM attachments.

2-storey top frame SM / TH: 2-storey top frame for all SM shakers. Dimensions adapted to TH 30. For vessels with a max. height of 14 cm (Erlenmeyer flasks 200 / 250 ml). Dimensions (w x d x h): 600 x 402 x 192 mm. For fastening the vessels, either universal trays SM or sliding plates with Combifix SM can be used.



Code	Description	Packaging
LLG06510572	Accessories for universal shaker SM 30, Type 2-storey top fram SM/TH (for using SM with TH)	
LLG09838211	Accessories for universal shaker SM 30, Type 2-storey top frame SM (without platforms)	1 pz.
LLG09838470	Incubator hood TH 30, Type Incubator Hood TH 30	1 pz.
LLG06255238	Accessories for universal shaker SM 30, Type Mounting set for universal trays or sliding plates	1 pz.
LLG09838201	Accessories for universal shaker SM 30, Type Rack system Combifix SM, Set A	1 pz.
LLG09838205	Accessories for universal shaker SM 30, Type Rack system Combifix SM, Set B	1 pz.
LLG09838208	Accessories for universal shaker SM 30, Type Rack system Combifix SM, Set C	1 pz.
LLG06231576	Accessories for universal shaker SM 30, Type Sliding plate for top frame with Combifix SM, Set A	1 pz.
LLG06285190	Accessories for universal shaker SM 30, Type Sliding plate for top frame with Combifix SM, Set B	1 pz.
LLG06285191	Accessories for universal shaker SM 30, Type Sliding plate for top frame with Combifix SM, Set C	1 pz.
LLG09838210	Accessories for universal shaker SM 30, Type Universal tray SM (without spring clamps)	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

