

**Combi shaker KL 2****Bühler****Combi shaker KL 2**

Mini shaker which can be switched between orbital and reciprocating motions. Includes flask attachment system with 3 clamping rollers. Versatile, for use in medical diagnostics, haematology, for cell cultivation or analysis. Max. load 4 x 500 ml or 2 x 1000 ml Erlenmeyer flasks.

Scope of supply: Shaker and accessory platform with 3 clamping rollers and anti-slip mat.

Speed range: 30 ... 420 rpm with stepless electronic adjustment

Orbit or stroke: 8 mm

Max. loading: 2 kg

Platform (W x D): 300 x 300 mm

Timer: 0 ... 120 min / continuous

Overall (W x D x H): 355 x 455 x 195 mm

Weight: 13,5 kg

Supply requirements: 230 V, 50/60 Hz

Protection class: IP 21



Code	Description	Packaging
LLG09838320	Combination shaker KL 2,reciprocat.motion speed: 30 - 420 rpm, platform:300 x 300 mm with attachment system with 3 clamping roll	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

