

**Incubator, mechanical control, B 28 series****BINDER****Incubator, mechanical control, B 28 series**

Robust, small, space-saving unit with mechanical controller and adjustable ventilation for accurate and reliable incubation.

- Temperature range: 30°C to 70°C
- Inner glass door
- Available with, or without, TB (Class 1)



Code	Description	Packaging
LLG09883556	Incubator B 28 without overheat cut-out, without TB, 28l, 230V, 50/60Hz	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

