

**Safety drying ovens, FDL****BINDER****Safety drying ovens, FDL**

A safety drying oven from BINDER ensures the perfect test result of your solvent-based paints and coatings, absolute temperature precision even with high ventilation.

- Defined solvent quantity according to EN 1539
- Electronically controlled APT.lineă preheating chamber
- Temperature range: from 5°C above ambient temperature to 300°C
- Independent adjustable temperature safety device class 2 (DIN 12880), with visual temperature alarm
- RS 422 interface for use with APT-COMă DataControlSystem communication software
- Power supply 230 V, 50/60 Hz

The number of shelves included in the scope of delivery is listed in the table.



| Code | Description | Packaging |
|-------------|--|-----------|
| LLG09883580 | Drying oven/lacquer drying oven FDL 115, 230V, 50/60Hz | 1 pz. |



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

