



## Portable autoclaves without heating

### Portable autoclaves without heating

For pressurised steam sterilisation. Up to 3 bar, 143°C. 4 bar version available on request. Chamber and lid made of stainless steel. Temperature indication by enclosed-scale thermometer.



Code	Description	Packaging
LLG09842061	Portable autoclaves without heating, Type LaS-3-20, Ext. diam. 400 mm, Height 590 mm, Int. diam. 260 mm, Int. height 400 mm	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

