



Graduated pipettes, piston type, AR-glas, white graduations

BRAND

**Graduated pipettes, piston type, AR-glas®, white graduations**

Silberbrand. Calibrated to deliver (TD, Ex), total delivery. Integral piston pipetting aid with plunger retaining spring.





**CARLO ERBA Reagents S.A.S.**  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

