

**Heating ovens, ED, FD, FED series - Classic.Line****BINDER****Heating ovens, ED, FD, FED series - Classic.Line**

Diversity for all types of thermal, whether efficient drying, long-term controlled elevated temperatures or sterilization tasks for homogeneous temperature distribution: a BINDER oven and heating chamber is up to any tasks thanks to its wide temperature range.

ED Series: Drying ovens with gravity convection.


FD Series: Drying ovens with forced convection.

FED Series: Heating chambers with forced convection and enhanced timer functions.

Equipment:

- Temperature range from 5°C above ambient temperature to 300°C
- Independent adjustable temperature safety device class 2 (DIN 12880), with visual temperature alarm
- Optional RS 422 Interface for APT-COMä Data Control System communication software

The number of shelves included in the scope of delivery is listed in the table.



Code	Description	Packaging
LLG09883803	Drying Ovens, Model FD400, upto +300°C, cap 400 ltrs	1 pz.
LLG09883529	General purpose ovens/hot Air sterilizers ED23, without RS 422, 23l	1 pz.
LLG09883612	Drying oven/hot-air sterilizer ED 23-UL, without RS 422, 23 ltr. with natural air-circulation, 115 V 1N ~ 60 Hz	
LLG09883615	Drying ovens Model FD 53, upto +300°C, cap. 53 ltrs	
LLG09883551	General purpose oven/hot air sterilizers with RS 422, ED 400	1 pz.
LLG09883808	Drying ovens FD 23, upto +300°C, 20 ltrs	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

