



Accessories for Spectrophotometer LLG-uniSPEC 2 / LLG-uniSPEC 4

Accessories for Spectrophotometer LLG-uniSPEC 2 / LLG-uniSPEC 4



Code	Description	Packaging
LLG06263617	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Cuvette holder, adjustable for 1 micro-cuvette	1 pz.
LLG06263622	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Cuvette holder, automatic for 8 cuvettes 10 mm	
LLG06266530	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Cuvette holder for 4 cuvettes 50 mm and 100 mm	1 pz.
LLG06263612	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Cuvette holder, water-cooled for 1 cuvette 10 mm	
LLG06263616	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Deuterium lamp	1 pz.
LLG06263623	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Peltier system	
LLG06266532	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Sipper system	
LLG06263615	Accessories for Spectrophotometer uniSPEC 2 / uniSPEC 4, Tungsten lamp	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

