



## Panoramic Eyeshield 6x3

### Panoramic Eyeshield 6x3

- CE approved to EN166 and EN170
- Clear polycarbonate lens with exclusive VANGUARD UDC coating (protection against scratching and fogging)
- Panoramic lens can be perfectly worn over prescription glasses
- Facial support compatible with half masks
- Additional facial protection available
- Indirect ventilation system assures certified protection against drops and splashes
- Replacement lens available
- Impact resistance at 432 km/h, even at extreme temperatures
- 100% UV protection (up to 380nm)
- Optimal wearing comfort thanks to the extra light weight of 75 g



Code	Description	Packaging
LLG06269305	Goggle PSA Cat. III UDC coating, indi/1	1 pz.
LLG06269306	Panoramic Eyeshield 6x3, Additional facial protection, Colour, Lens clear, polycarbonate, 2C-1.2	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

