

**pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20**

Mettler-Toledo Online

pH meter FiveEasyä pH F20 / FiveEasy Plusä pH FP20

FiveEasyä and FiveEasyä Plus benchtop meters provide high quality pH/mV measurements with the simple click of a button.

- Large, easy to read and easy to understand display with all information visible at a glance
- Five self-explanatory buttons making operation simple and easy
- Fast and easy connectivity for simple and secure data transfer (FP only)
- Electrode arm extension pole included

Scope of supply:

F20-Meter: Meter with CD operating instructions, QuickGuide, declaration of conformity, test certificate, power supply

F20-Standard: as F20-Meter, plus plastic LE438 3-in-1 pH electrode, 1 each x 4.01 pH, 7.00 pH, 9.21 pH buffer sachets

FP20-Meter: Meter with CD operating instructions, QuickGuide, declaration of conformity, test certificate, power supply

FP20-Standard: as FP20-Meter, plus plastic LE438 3-in-1 pH electrode, 1 each x 4.01 pH, 7.00 pH, 9.21 pH buffer sachets

FP20-Bio: as FP20-Meter, plus glass LE410 3-in-1 pH electrode, 1 each x 4.01 pH, 7.00 pH, 9.21 pH buffer sachets

FP20-Micro: as FP20-Meter, plus glass LE422 micro pH electrode, 1 each x 4.01 pH, 7.00 pH, 9.21 pH buffer sachets

FP20-TRIS: as FP20-Meter, plus glass LE420 pH electrode with moveable sleeve diaphragm, 1 each x 4.01 pH, 7.00 pH, 9.21 pH buffer sachets

Specifications**Measuring range (Accuracy)**

pH F20 / FP20: 0.00 ... 14.00pH / -2.00 ... 16.00pH (±0.



01pH)

mV: -2000 ... 2000mV (± 1 mV)

Temperature F20 / FP20: 0 ... 100°C (0.5°C) / -5 ... 105°C (0.3°C)

Input: BNC, Cinch (NTC 30kW)

Temperature compensation: ATC and MTC

Power supply: 100 ... 240V / 50 ... 60Hz / 12V DC

Calibration F20: 3 points, linear, 4 predefined buffer groups

Calibration FP20: 5 points, linear/segmented, 4 predefined buffer groups

FP20 only

Memory: 200 measurements

Interface: RS232, USB, Reference input

Code	Description	Packaging
LLG06280100	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type F20-Meter	1 pezzo
LLG06280101	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type F20-Standard Kit	1 pz
LLG06280104	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type FP20-Bio	1 pz.
LLG06280102	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type FP20-Meter	1 pezzo



Code	Description	Packaging
LLG06280105	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type FP20-Micro	1 pz.
LLG06280103	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type FP20-Standard Kit	1 pezzo
LLG06280106	pH meter FiveEasy pH F20 / FiveEasy Plus pH FP20, Type FP20-TRIS	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

