

**Conductivity meters FiveEasy F30 / FiveEasy Plus FP30**

Mettler-Toledo Online

**Conductivity meters FiveEasyä F30 / FiveEasyä Plus FP30**

FiveEasyä and FiveEasyä Plus benchtop meters provide high quality conductivity measurements with the simple click of a button. They are the perfect choice for reliable measurement data at a reasonable price.

- Large, easy to read and easy to understand display with all information visible at a glance
- Five self-explanatory buttons making operation simple and easy
- Fast and easy connectivity for simple and secure data transfer (FP only)
- Electrode arm extension pole included

**Scope of supply:**

**F30/FP30-Meter:** Meter with CD operating instructions, QuickGuide, declaration of conformity, test certificate, power supply

**F30/FP30-Standard:** as F30/FP30-Meter, plus LE703 conductivity sensor, 1 each x 1413  $\mu\text{S}/\text{cm}$  and 12.88  $\text{mS}/\text{cm}$  standard sachets

**Specifications****Measuring range (Accuracy)**

**Conductivity F30 / FP30:** 0.01  $\mu\text{S}/\text{cm}$  ... 200  $\text{mS}/\text{cm}$  / 0.01  $\mu\text{S}/\text{cm}$  ... 500  $\text{mS}/\text{cm}$  ( $\pm 0.5\%$  of measured value)

**TDS F30 / FP30:** 0.01  $\text{mg}/\text{l}$  ... 200  $\text{g}/\text{l}$  / 0.01  $\text{mg}/\text{l}$  ... 300  $\text{g}/\text{l}$

**Temperature F30 / FP30:** 0 ... 100  $^{\circ}\text{C}$  (0.5  $^{\circ}\text{C}$ ) / -5 ... 105  $^{\circ}\text{C}$  (0.3  $^{\circ}\text{C}$ )

**Calibration F30:** 1 point, 3 predefined standards

**Calibration FP30:** 1 point, 3 predefined standards, manual cell constant input

**Input:** Mini-DIN

**Temperature compensation:** Linear: 0.00  $\%/^{\circ}\text{C}$  ... 10.00  $\%/^{\circ}\text{C}$

**Reference temperature:** 20  $^{\circ}\text{C}$  and 25  $^{\circ}\text{C}$

**Power supply:** 100 ... 240 V / 50 ... 60 Hz / 12 V DC

**FP30 only**

**Salinity:** 0.00 ... 42 psu



**Memory: 200 measurements**

**Interface: RS232, USB**

Code	Description	Packaging
LLG06280107	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type F30-Meter	1 pz.
LLG06280108	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type F30-Standard Kit	1 pz
LLG06280109	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type FP30-Meter	1 pz.
LLG06280110	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type FP30-Standard Kit	1 pz





CARLO ERBA Reagents S.A.S.  
Chaussée du Vexin – 27106 Val  
de Reuil cedex  
N° TVA: FR 63391048824  
Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



*CARLO ERBA  
Reagents operates with  
a Certified Quality  
Management System*

