

**Conductivity meters FiveEasy F30 / FiveEasy Plus FP30**

Mettler-Toledo Online

Conductivity meters FiveEasyä F30 / FiveEasyä Plus FP30

FiveEasyä and FiveEasyä Plus benchtop meters provide high quality conductivity measurements with the simple click of a button. They are the perfect choice for reliable measurement data at a reasonable price.

- Large, easy to read and easy to understand display with all information visible at a glance
- Five self-explanatory buttons making operation simple and easy
- Fast and easy connectivity for simple and secure data transfer (FP only)
- Electrode arm extension pole included

Scope of supply:

F30/FP30-Meter: Meter with CD operating instructions, QuickGuide, declaration of conformity, test certificate, power supply

F30/FP30-Standard: as F30/FP30-Meter, plus LE703 conductivity sensor, 1 each x 1413 $\mu\text{S}/\text{cm}$ and 12.88 mS/cm standard sachets

Specifications**Measuring range (Accuracy)**

Conductivity F30 / FP30: 0.01 $\mu\text{S}/\text{cm}$... 200 mS/cm / 0.01 $\mu\text{S}/\text{cm}$... 500 mS/cm ($\pm 0.5\%$ of measured value)

TDS F30 / FP30: 0.01 mg/l ... 200 g/l / 0.01 mg/l ... 300 g/l

Temperature F30 / FP30: 0 ... 100 $^{\circ}\text{C}$ (0.5 $^{\circ}\text{C}$) / -5 ... 105 $^{\circ}\text{C}$ (0.3 $^{\circ}\text{C}$)

Calibration F30: 1 point, 3 predefined standards

Calibration FP30: 1 point, 3 predefined standards, manual cell constant input

Input: Mini-DIN

Temperature compensation: Linear: 0.00 $\%/^{\circ}\text{C}$... 10.00 $\%/^{\circ}\text{C}$

Reference temperature: 20 $^{\circ}\text{C}$ and 25 $^{\circ}\text{C}$

Power supply: 100 ... 240 V / 50 ... 60 Hz / 12 V DC

FP30 only

Salinity: 0.00 ... 42 psu



Memory: 200 measurements

Interface: RS232, USB

Code	Description	Packaging
LLG06280107	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type F30-Meter	1 pz.
LLG06280108	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type F30-Standard Kit	1 pz
LLG06280109	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type FP30-Meter	1 pz.
LLG06280110	Conductivity meters FiveEasy F30 / FiveEasy Plus FP30, Type FP30-Standard Kit	1 pz





CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

