

**LLG-Disposable Protective Visors****LLG Labware****LLG-Disposable Protective Visors**

- Transparent visor made of acetate
- With anti-fog coating
- Head holder made of TPU-foam with elastic
- Ideal as a disposable protection, for example, for visitors
- Can be worn over prescription glasses and protective masks
- Weight: 22 g



Code	Description	Packaging
LLG06284812	Disposable Protective Visors, Disposable Protective Visors, clear lens, elastic headband, anti-fog	20 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

