

**Metal trays for Histology****Kartell****Metal trays for Histology**

Made of stainless steel for easy removal of paraffin inclusions: they allow a precise location of the specimen into the paraffin.



Code	Description	Packaging
LLG06285348	Metal trays for Histology, External dim. (W x D x H) 52 x 35 x 11 mm	10 pz.
LLG06285347	Metal trays for Histology, External dim. (W x D x H) 52 x 35 x 11,6 mm	10 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

