



Rotors for Concentrator plus

Rotors for Concentrator plusä



Code	Description	Packaging
LLG09776493	Rotors for Concentrator plus, Type F-40-18-19, for 18 tubes 10 ml (19 x 66 mm)	
LLG09776489	Rotors for Concentrator plus, Type F-45-24-12, for 24 tubes 5 ... 7 ml (12 x 67 ... 100 mm)	
LLG09776491	Rotors for Concentrator plus, Type F-40-36-12, for 36 tubes 1.5 ml (12 x 32 mm)	
LLG09776492	Rotors for Concentrator plus, Type F-45-36-15, for 36 tubes 3.0 / 5.0 ml (15 x 45 ... 48 mm)	
LLG09776465	Rotors for Concentrator plus, Type F-45-48-11, for 48 tubes 1.5 / 2.0 ml (11 x 41/11 x 47 mm)	1 pz.
LLG09776499	Rotors for Concentrator plus, Type F-45-70-11, for 70 tubes 1.5 / 2.0 ml (11 x 41/11 x 47 mm)	
LLG09776498	Rotors for Concentrator plus, Type F-45-72-8, for 72 tubes 0.5 ml (8 x 31 mm)	
LLG09776463	Rotors for Concentrator plus, Type F-50-8-16, for 8 tubes 8 ... 12 ml (16 x 105 ... 120 mm)	
LLG09776496	Rotors for Concentrator plus, Type A-2-VC, for MTP, PCR and Deepwell plates (128 x 86 x 27 mm)	1 pz.
LLG09776624	Rotor Type F-35-6-30 , for 6 conical tubes 15 / 50 ml (17 / 30 x 116 ... 123 mm)	1 pz.
LLG09776494	Rotors for Concentrator plus, Type F-35-8-24, for 8 tubes 25 ml (24 x 86 ... 105 mm)	
LLG09776490	Rotors for Concentrator plus, Type F-45-8-17, for 8 conical tubes 15 ml (17 x 116 ... 123 mm)	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

