



LED Illumination, direct

SCHOTT

LED Illumination, direct

EasyLED Ringlight / EasyLED Ringlight Plus: Powerful and robust LED ring light with continuous dimming. Ringlight Plus version with basic segment features and segment rotation. Easy and intuitive operation via integrated jog dial wheel. Including wide input range power supply with world wide plug (clip system). External / internal diameter 114/66 mm. If necessary, please order the necessary lens adapter separately.

EasyLED Spotlight Plus (Single/ Double): Advanced high power LED gooseneck system with integrated controller and continuous dimming. With the Spotlight Plus Double, the two spots can also be activated individually. Including wide input range power supply with world wide plug (clip system). **Please order mounting bracket or column adapter separately** (same adapter as for KL 300 LED).



Code	Description	Packaging
LLG06238244	LED Illumination, direct, Type Ringlight Plus, 8 segments, Working Distance 65 mm, Colour temperature 5600 K, Light flux 250 lm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

