

**Antibiotic test papers****Hahnemühle****Antibiotic test papers**

Recommended for identification of pathogens of infectious diseases by determination of the degree of resistance against antibiotics

- Made from ultrapure raw materials, entirely without additives to ensure that no inhibition will occur during the incubation
- Consistent thickness and absorption volume
- Available as small circles, even individually printed on one or both sides
- Meet the directive for carriers for active agents acc. to DIN 58940-2
- Further sizes on request



Code	Description	Packaging
LLG04006131	Antibiotic test papers, Type 3324, Diam. 6 mm, Thickness 0.73 mm, Weight 280 g / m ²	1000 pz.
LLG04006126	Antibiotic test papers, Type 22, Diam. 6 mm, Thickness 0.35 mm, Weight 180 g / m ²	1000 pz.
LLG04006129	Antibiotic test papers, Type 2668, Diam. 6 mm, Thickness 0.9 mm, Weight 320 g / m ²	1000 pz.
LLG04006127	Antibiotic test papers, Type 22, Diam. 9 mm, Thickness 0.35 mm, Weight 180 g / m ²	1000 pz.
LLG04006130	Antibiotic test papers, Type 2668, Diam. 9 mm, Thickness 0.9 mm, Weight 320 g / m ²	1000 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

